

CENTRIVAP™ CONSOLE

78120-00 (115 V)
78120-01 (230 V)

INSTRUCTION MANUAL

Product designs are subject to change without notice
© 1993 Labconco Corporation
Form 78303 REV A / ECO 7320
Printed in U.S.A.

TABLE OF CONTENTS

Introduction

Components Shipped	3
General Description	4
Component Identification	6 - 9

Installation

Electrical Connections	10
Vacuum Pump	10
Rotor Installation	11
Cold Trap Installation	12
Secondary Trap	13

Normal Operation

CentriVap Console	14
Auxiliary Equipment	15

Applications

16

Routine Maintenance Schedule

17

Specifications

18

Replacement Parts

20-23

Dimensions

24

Schematics

25-26

Wiring Diagram

27-28

Vacuum Plumbing Diagram

29

Troubleshooting

30-31

Warranty

32

Shipping Claims

33

Accessories

34-35

Contacting Labconco

36

Components Shipped

Carefully check the contents of the carton for damage that might have occurred in transit. Do not discard the carton or packaging material until all components have been checked against the following component list and the equipment has been installed and tested.

As shipped, the carton should contain the following:

Part Number	Description
78120-00 or 78120-01	CentriVap Console (115 V/60 Hz) CentriVap Console (230 V/50 Hz)
78303	Instruction Manual
78164	Rotor for (54) 1.5 ml Vials
78152	Solvent Trap Insert

General Description

The CentriVap Console combine centrifugal force with vacuum to rapidly evaporate solvents from biological or analytical samples. Centrifugation eliminates bumping and foaming as vacuum is applied and concentrates the solute in the bottom of the vial. This allows recovery of solutes from sample volumes as small as a few microliters. The unit is equipped with a 300 watt heater to speed evaporation by warming the chamber to 40°C during processing. An insulating fiberglass barrier reduces heat generated by magnetic coupling from being transferred into the sample chamber. When the heat mode is selected, the start up of the vacuum pump is delayed 5 minutes to allow more efficient heat transfer to the samples. A hexagonal rotor capable of holding (54) 1.5 ml vials is included. A second rotor capable of holding (18) 1.5 ml vials may be added. Optional rotors for 12mm x 75mm tubes and 50ml centrifugal tubes are also available.

The built-in cold trap protects the vacuum pump by trapping moisture, vapors and corrosive fumes as they evaporate from the samples. The high speed condensing unit and the low temperature refrigerant reduce the temperature of the cold trap to -60°C in 30 minutes. The stainless steel trap is used for aqueous and organic solvents. The optional glass trap insert should be used for corrosive solvents. A secondary trap is supplied to protect the vacuum pump (not supplied) from most vapors which might pass through the cold trap. The CentriVap is supplied with a cartridge for trapping solvents. Other cartridges are available for trapping water, acids and radioisotopes.




Figure 1

Component Identification

- (1) **Refrigeration Switch.** Activates the refrigeration system that cools the cold trap.
- (2) **Rotor Switch.** Activates the rotation of the rotor and turns on the vacuum pump.
- (3) **Heater Switch.** Activates the heater used for warming the samples. When both the rotor and heater switches are turned on, the start up of the vacuum pump is delayed 5 minutes to allow time for the samples to be efficiently warmed.
- (4) **Vacuum Switch.** Activates the vacuum pump independently from the rotor operation control.
- (5) **Vacuum Gauge.** Indicates the system vacuum in microns when the vacuum pump is operating.
- (6) **Vacuum Control Knob (labeled 'Concentrator')** - Rotating the vacuum control knob so the indicator points to "Vacuum" directs vacuum to the CentriVap chamber. Rotating the vacuum control knob so the indicator points to slightly below horizontal isolates the chamber from the vacuum pump and maintains the existing vacuum in the chamber. Rotating the vacuum control knob so the indicator points to "Vent/Backfill" releases the vacuum in the chamber.
- (7) **Vent/Backfill Port.** The chamber vacuum is released through this port. An inert gas may be introduced into the chamber through this port when the Vacuum Control Knob is positioned to the Vent/Backfill position.
- (8) **Vacuum Control Knob (labeled 'Auxiliary')** - Rotating the vacuum control knob so the indicator points to "Vacuum" directs vacuum to the auxiliary vacuum port on the rear of the cabinet. Rotating the vacuum control knob so the indicator points to slightly below horizontal isolates the auxiliary vacuum port from the vacuum pump and maintains existing vacuum in any device connected to the auxiliary vacuum port. Rotating the vacuum control knob so the indicator points to "Vent" releases vacuum in any device that is connected to the auxiliary vacuum port.
- (9) **Clear Acrylic Lid.** Permits samples to be examined without opening the lid. Never attempt to open the lid unless the rotor is completely stopped.
- (10) **Rotor.** Holds sample vials or tubes during concentration. The rotor is belt driven for smooth, quiet operation and is magnetically coupled.

Component Identification (Cont)

- (11) **Stainless Steel Trap.** Cooled to -60°C , this trap collects the solvent that is evaporated from the sample. When evaporating corrosive chemicals the optional glass trap should be inserted into the stainless steel trap.
- (12) **Secondary Trap.** Protects the vacuum pump by trapping most vapors which may pass through the Cold Trap. The CentriVap is equipped with a cartridge for trapping solvents. Other cartridges are available for water, acids and radioisotopes.
- (13) **Vacuum Pump Receptacle.** Provides electrical power for the vacuum pump. Pump operation is delayed 5 minutes when the heater mode is selected.
- (14) **Auxiliary Vacuum Port.** Provides a means to connect auxiliary equipment to the primary cold trap, secondary trap and vacuum source. The Vacuum Control Knob (labeled 'Auxiliary') controls vacuum to this port. To capture vapors in the primary cold trap, the refrigeration system must be activated.
- (15) **Quick Disconnect Fittings.** Provides a quick and easy means to disconnect hoses from the glass trap accessory at the end of a run.
- (16) **Neoprene Gasket.** Provides a sealing surface between the chamber and acrylic top.
- (17) **300 W Heater.** Provides heat to chamber controlled at 40°C when activated.
- (18) **Thermally Protected Motor.** Provides protection from motor overheating.
- (19) **Lid Latch.** A latch prevents the lid from being opened while the rotor is turning.
- (20) **Belt Drive System.** Rotates rotor up to 1800 RPM. Ensures smooth, quiet operation and consistent, even braking action.
- (21) **Optional Glass Trap Insert.** For corrosive systems, this optional glass insert is used within the stainless steel trap.
- (22) **Condensing Unit.** High speed 1/4 HP refrigeration system paired with the low temperature refrigerant reduces temperature of the Cold Trap to -60°C .

Component Identification (Con't)

- (23) **Acrylic Lid - Cold Trap.** Used specifically with the stainless steel trap. Fittings are supplied for 1/2" internal diameter (ID) tubing.
- (24) **Cast Aluminum Chamber.** Double impregnated for vacuum integrity and coated for easy cleaning of any spills. Large volume accommodates single or double deck rotors.

A. Electrical Connections

The Labconco CentriVap Console, Model 78120-00, is fully wired with an electrical cord and plug to connect to a standard outlet. For this model, a fused 115 V, 60 Hz, 1-phase, 20 amp supply is required. The Labconco CentriVap Console, Model 78120-01, is wired without a plug on the end of the electrical cord. For this model, a fused 230 V, 50 Hz, 1-phase, 10 amp supply and plug are required.

B. Vacuum Pump

A vacuum pump must be supplied by the user. Generally, the lower the vacuum that the pump can produce and the greater the pumping capacity, the better the performance will be that can be expected from the CentriVap Console. In order to achieve very low vacuums, usually rotary vane direct drive vacuum pumps are used. If oil contamination of the rotary vane pump is considered objectional, an alternative is the use of a corrosion resistant diaphragm pump capable of pumping at least 74 liters/minute and obtaining at least 29 inches of mercury vacuum. See Accessories Pages 34 thru 35 for pump selections. To install the vacuum pump in the CentriVap console:

- (1) Remove the front access panel by unscrewing the two black fasteners located in the upper corners.
- (2) The vacuum pump fitting must be capable of accepting 1/2 inch ID tubing.
- (3) Position the vacuum pump on the lower shelf so that the vacuum port is beneath the clearance hole in the right side of the middle shelf.
- (4) Feed the existing loose tube coming from the secondary trap through the hole in the shelf and attach it to the vacuum port of the pump. Install the hose clamp.
- (5) Plug the vacuum pump line cord into the receptacle at the back of the cabinet. If the pump is equipped with a power switch, then turn the pump power switch to "ON".

C. Rotor Installation

- (1) Unpack and clean the rotor(s). If the dual rotors are used, install the smaller rotor on the shaft and rotate until the slots on the bottom of the rotor engage the pin in the shaft.
- (2) Place the top rotor over the shaft. Gently push down on the rotor button once to engage and hold the rotor (Figure 2)
- (3) If the rotor does not have the chrome colored button in the top, place the rotor over the shaft and lower it into position. Rotate it until the slots on the bottom engage the pin in the shaft.




Figure 2

D. Cold Trap

- (1) If your application involves corrosive chemicals, go to step 4, skipping steps 2 and 3.

- (2) Place the acrylic top fitted with the gasket insulation and quick disconnects over the stainless steel trap.

- (3) Turn Cold Trap switch to "On" position. Allow at least 30 minutes for the cold trap to come to temperature. Caution: Labconco recommends the use of the optional glass trap insert for most applications. Please read the following section even if your work does not involve corrosive chemicals.

- (4) If using corrosive chemicals and the optional glass trap, measure approximately 200ml ethanol (ethyl alcohol) as a heat transfer medium into the stainless trap.

- (5) Insert the glass trap into the ethanol filled trap and check to see if the ethanol level is at least 2/3 up the outside of the glass. If not, add ethanol until this level is reached.

- (6) Fit the glass trap insulation over the two ports in the glass trap and press the insulation down into the console top.

- (7) Install 90° elbows onto the two glass ports by slowly and gently pushing them on.

CAUTION: Excessive force could break the glass.

- (8) Remove the two hoses and barb fittings from the existing acrylic cold trap cover by pressing the dark gray collar in and simultaneously pulling the barb fitting out.

- (9) Push the two barb fittings into the elbows on the top of the glass trap. The left end of the front hose should connect to the center elbow of the glass trap.

- (10) Turn Cold Trap switch to "On" position. Allow at least 30 minutes for the cold trap to come to temperature.

- (11) Immediately following a run with the optional glass trap, defrost by running under cold water.

CAUTION: Failure to immediately defrost could break the glass trap.

E. Secondary Trap

The Labconco CentriVap Console comes with a secondary trap installed in place and a solvent trap insert which must be installed. If acidic, aqueous or radioactive samples are to be evaporated, order the appropriate insert. Then install the chosen insert as follows:

- (1) Unscrew the clear bottom of the trap.
- (2) Pull any existing cartridge down from the head to remove it.
- (3) Press the neck of the new insert into the head until both O-rings are inserted.
- (4) Reinstall the clear bottom of the trap by screwing it tightly to the head.
- (5) Properly dispose of used inserts as governed by federal, state or local regulation.

After a new insert is installed it may be necessary to run the vacuum pump for an extended period of time before the reading on the vacuum gauge will reproduce previous low readings. This procedure causes moisture to be purged from the new insert. Once this procedure is completed, subsequent runs should operate in the normal manner.

CentriVap Console

- (1) Turn the CentriVap Cold Trap switch to "ON" for 30 minutes to cool down to operating temperatures. Do not turn on the vacuum or rotor at this time.
- (2) Load samples into rotors(s).
- CAUTION: Do not operate with an unbalanced rotor.**
- (3) Close lid.
- (4) The vacuum pump should be connected to the timed outlet on the back of the CentriVap Console if the delay feature is desired.
- (5) Rotate the Concentrator Vacuum Control Knob to "VACUUM" and the Auxiliary Vacuum Control Knob to "VENT".
- (6) If heat is desired, turn the Heater switch to the "ON" position before turning on the rotor switch. If heat is chosen, the vacuum pump will be delayed for 5 minutes to provide adequate heating to the samples prior to evacuation. If heat is not chosen, the vacuum pump will start when the rotor switch is turned on.
- (7) Turn rotor switch to "ON". If heat was not selected, the vacuum pump will start immediately. If heat was selected, the vacuum pump will start after a five minute delay.
- (8) At the end of a run, vent the vacuum by turning the Concentrator Vacuum Control Knob to "VENT" and turn the rotor switch to "OFF". Do not lift the lid until the rotor is completely stopped and the vacuum is vented.
- (9) Immediately following a run using the optional glass trap, defrost by running the glass trap under cold water.
- CAUTION: Failure to immediately defrost could break the glass trap.**
- (10) At the end of the day, the vacuum line to the pump may also be vented by either turning the Concentrator Vacuum Knob to "VACUUM" while the lid is open or by turning the Auxiliary Vacuum Knob to "VACUUM" when nothing is attached to the auxiliary port at the rear of the cabinet.

Auxiliary Equipment

- (1) Attach the auxiliary equipment to the auxiliary vacuum port on the rear of the cabinet with a vacuum hose securing the vacuum hose with hose clamps.
- (2) Turn the CentriVap Cold Trap switch to "On" for 30 minutes to cool down to operating temperature. Do not turn on the vacuum at this time.
- (3) Rotate the Concentrator Vacuum Control Knob to "Vent" and the Auxiliary Vacuum Control Knob to "Vacuum".
- (4) After the CentriVap Cold Trap switch has been "On" for 30 minutes, the CentriVap Vacuum switch may be turned "On".
- (5) The auxiliary equipment may now be run.
- (6) At the end of an auxiliary vacuum run, vent the vacuum by rotating the Auxiliary Vacuum Control Knob to "Vent".

- A. Hazardous Solvents/Ammonium Hydroxide.** Most of these will condense in the Cold Trap. Install a secondary trap insert for organic solvents (P/N 78152) to capture any vapors and protect the pump. Handle these tubes with care.

- B. Radioactive Substances**
When evaporating radioactive substances, install an activated carbon insert (P/N 78150) in the secondary trap and use in conjunction with the Cold Trap. The carbon insert is intended to protect the vacuum pump from traces of radioactive materials that may have passed through the cold trap. It does not meet NRC filter design recommendations. After operating, properly dispose of the insert as radioactive waste.

- C. Volatile Organics**
Organic solvents, such as acetone and chloroform, are very volatile. The Cold Trap operating at -60°C might not be sufficient to collect 100 percent of the vapor; therefore, use of a diaphragm pump which will produce $29''\text{ Hg}$ vacuum to evaporate these solvents is required.

CAUTION: When corrosive or volatile organics are used, a secondary trap insert must be placed between the diaphragm vacuum pump and Cold Trap to protect the pump.

- D. Corrosive Solvents**
As a precaution, install an acid insert (P/N 78148) in the secondary trap canister to minimize any acid vapors from entering the pump.
- E. Pump Exhaust Filter**
Install a pump inlet filter (P/N 14723) on the pump to eliminate any mists.

- F. Clean Room Applications**
To ensure all exhaust from the pump is eliminated from the room, attach tubing to the oil mist eliminator and route to a vented hood.

- G. Routine Aqueous Applications**
Install a secondary trap insert for water vapor (P/N 78149) to protect the pump.
- NOTE:** Dispose of all chemicals and hazardous materials in compliance with local codes.

Labconco is not responsible for illegal dumping of any chemicals or materials.

Concentrator And Cold Trap

Before any maintenance is performed, all power switches should be "Off" and the CentriVap should be unplugged from electrical outlets.

- (1) Do not store any material around or on top of the CentriVap unit. Proper air ventilation is essential.
- (2) Clean any spills immediately to prevent solvents from accumulating on rotors and inside the chamber.
- (3) Check the oil level of the pump if applicable. It should be between "MIN" and "MAX". If the oil level is low, add oil to proper level.
- (4) If the vacuum pump oil shows cloudiness, particles or discoloration, drain the pump and replace with fresh oil.
- (5) Check the Cold Trap for condensed or frozen solvents and dispose of them appropriately. Completely empty the trap before the next run. If solvents are frozen in the glass trap, run it under cold water immediately after operating.
- (6) If the glass trap is used, check to see that the ethanol in the stainless steel trap is free of ice water. Aspirate or pipet the contaminated ethanol and replace it with 200 ml of fresh ethanol.
- (7) If the media in the secondary trap has completely changed color, discard and replace with a new insert. For the activated carbon insert used in radioactive applications, no indicator exists. In radioactive applications, the system should be monitored with a Geiger counter.

NOTE: The user is solely responsible for any radioactive leak.

- (8) Recheck all connections for a secure fit to prevent leakage.

A. CONCENTRATOR

Heater: 300W, built-in tubular heater controlled at 40°C with over-temperature protection
Power Input: 115 VAC, 1 phase, 60 Hz, 8.5 amp (Model 78120-00) 230 VAC, 1 phase, 50 Hz, 4.3 amp (Model 78120-01)
Drive: Belt drive coupled with magnetic drive
Controls: Rocker switches & rotary valves
Motor Rotation: Up to 1800 RPM
Lid: Abrasion resistant acrylic
Chamber: Double impregnated, epoxy coated aluminum casting
Cabinet: Epoxy coated steel
Dimensions: 36 3/4"H x 24"D x 18"W (93.3cm x 61.0cm x 45.7cm)
Shipping Weight: 167 lb. (75.7 kg)
Vacuum Gauge: Atmosphere to 5 microns
Rotors: (54) 1.5 ml vials (std) See Accessories Pages 34 -35 for other rotors.

B. PRIMARY COLD TRAP

Temperature: -60°C ±5°C @ 24°C ambient
Refrigeration: 1/4 HP condensing unit
Trap: 1 liter liquid, 500ml ice
Glass Trap Insert 375 ml liquid
(Optional)

C. SECONDARY TRAP

Inserts: Solvent (Std)
Water (Opt)
Acid (Opt)
Radioisotope (Opt)




Figure 3

REPLACEMENT PARTS

ITEM	REQ'D	MODEL 78120-00	MODEL 78120-01	DESCRIPTION
		PART NO.	PART NO.	
1	1	78255	78255	Thermostat Temperature Controller
2	1	78254	78254	Thermostat Temperature Limiter
3	1	78230	78230	Heater
4	1	78252	78252	Belt-Drive
5	1	78102	78102	Post Support Assembly
6	1	78250	78250	O-Ring
7	1	78524	78524	Lid
8	1	78297	78297	Motor Assembly
9	1	78295	78295	Capacitor - Motor
10	1	78525	78525	Gasket - Concentrator
11	1	78557	78557	Insulation - Trap
12	1	78290	78290	Timer - Delay
13	1	13498	N/A	Fuse 20A
14	2	N/A	13352	Fuse 10A
15	1	77343	77344	Condensing Unit
16	1	78125	78126	Condensing Module Ass'y (Inc. #15 & 17 & 44 & 46)
17	11 oz.	76222-01	76222-01	Refrigerant 13 B1
18	4	19414	19414	Caster
19	1	78129	78129	Secondary Trap Assembly
20	1	78506	78506	Lid Cold Trap
21	1	78505	78505	Gasket - Cold Trap
22	4	78286-01	78286-01	Hose Clear 24" Long
23	7	14888	14888	Hose Clamp
24	1	78370	78370	Speed Sensor Ass'y.
25	1	78148	78148	Acid Trap Insert
26	1	78149	78149	Water Vapor Trap Insert
27	1	78152	78152	Organic Solvent Trap Insert
28	1	78150	78150	Radioisotope Trap Insert
29	1	78156	78156	Glass Trap
30	1	78197	78197	Insulator
31	2	15445-01	15445-01	Elbow - Union
32	4	15481-01	15481-01	Stem - Tube to Hose
33	2	15487-01	15487-01	Bulkhead - Union
34	2	15489-01	15489-01	Elbow - Plug In
35	1	78164	78164	Rotor 54 x 1.5 ml Microcentrifuge
36	1	78169	78169	Rotor 18 x 1.5 ml Microcentrifuge (lower tier)
39	1	78512	78512-01	Solenoid
40	1	13019	13019	Switch - Lid
41	1	78528	78528	Switch, Magnetic Reed
42	1	12848	13339	Receptacle
43	1	78218	78218	Plate, Magnet Cover
44	1	78179	78179	Evaporator Coil
45	1	78292	78292	Cast Aluminum Chamber
46	1	14895	14895	Refrigeration Filter





Figure 4

REPLACEMENT PARTS

ITEM	REQ'D	PART NO.	DESCRIPTION
1	1	76815	Vacuum Gauge
2	2	13015-00	Switch - Cold Trap, Rotor
3	2	13015-01	Switch - Heater, Vacuum
4	2	75900	Valve Assembly
5	2	78333	Knob
6	4	78273	Lamp Assembly

Figure 5



Model 78120-00
(115 V)




Figure 6




Figure 7

WIRING DIAGRAM

Model 78120-00
(115V)




Figure 8

Model 78120-01
(230 V)





Figure 9

VACUUM PLUMBING DIAGRAM



PART NO. 78573

Figure 10

PROBLEM	CAUSES	CORRECTIVE ACTION
Lack of adequate vacuum. Gauge reads above 200 microns when chamber is empty (unloaded)	Air leak in the tubing	Check tubing for wear and replace. Tighten all connections
	Pump is not functioning properly	Check pump by locating a vacuum gauge closer to pump and close off rest of system. Check pump oil for cloudiness or particles and change. If pump is faulty, seek authorized service or replace pump
	Concentrator gasket is not sealing properly	Check gasket for cleanliness Adjust gasket as needed
	Secondary trap canister not fitted properly	Tighten all connections to and from the secondary trap. Tighten clear plastic housing.
	Secondary trap insert is spent	Replace with new insert
	Pump is not capable of producing deep vacuum	Change vacuum pump or continue operating with existing pump (Diaphragm pumps are not capable of producing deep vacuums)
Samples take too long to evaporate	Collapsed tubing	Check for collapsed tubing and replace
	Vacuum reading at concentrator with no load is greater than 200 microns	Check for seal of gasket and resolve or replace

PROBLEM	CAUSES	CORRECTIVE ACTION
Samples take too long to evaporate (Con't)	Obstruction in line	Trace by blocking various positions and reading gauge which is now placed right before vacuum pump. Remove obstructed tubing and replace
	Quick disconnects properly seated	Push quick disconnects to seal
Recovery of condensate in cold trap is less than normally expected	Cold trap is not ON	Check to make sure switch is ON and condensing unit fan is moving air out of the rear of the cabinet.
	Cold trap does not cool down	Turn vacuum pump OFF and allow Cold Trap to cool for at least 30 minutes to reach temperature -60° temperature can be checked with a solvent thermometer or digital thermometer.
Frequent oil change needed in pump	Secondary trap insert is spent	Change insert often
	Cold trap is not emptied after each run and dried	Empty the traps (glass or stainless steel) after each run and replace
	Vacuum too strong for chemical	Use a secondary trap insert and diaphragm pump to pull 29" Hg.

We are committed to providing our customers with quality equipment and service after the sale. Part of this objective involves keeping you informed of changes and new product additions. We therefore request that you take a moment to fill out the product registration card so we may know your location as well as some of the reasons that prompted you to purchase our products.

Labconco Corporation warrants products of its manufacture for one year, from receipt of the equipment by the purchaser, against defects in materials and workmanship. This limited warranty covers parts and labor but not transportation and insurance charges. In the event of a warranty claim contact the dealer who sold you the product. If the cause is determined to be a manufacturing fault, the dealer or Labconco Corporation will repair or replace all defective parts to restore the unit to operation. Under no circumstances shall Labconco Corporation be liable for indirect, consequential or special damages of any kind. This statement of warranty may be altered by a specific published amendment. No individual has authorization to alter the provisions of this warranty policy or its amendments. Lamps and expendable items such as filters are not covered by this warranty. Damage due to corrosion or accidental breakage are also not covered.

WARNING: The disposal and/or emission of substances used in connection with this equipment may be governed by various federal, state or local regulations. All users of this equipment are urged to become familiar with any regulations that apply in the user's area concerning the dumping of waste materials in or upon water, land or air and to comply with such regulations.

If a shipment is received in visibly damaged condition, be certain to make a notation on the delivering carrier's receipt and have his agent confirm the damage on your receipt. Otherwise, the damage claim may be refused.

If concealed damage or pilferage is discovered, notify the carrier immediately and retain the entire shipment intact for inspection. Interstate Commerce Commission rules require that the claim be filed with the carrier within 15 days after delivery.

NOTE: Do not return goods. Goods returned without prior authorization will not be accepted. Labconco Corporation and its dealers are not responsible for shipping damage. Claims must be filed directly with the freight carrier by the recipient. If authorization has been received to return this product, by accepting this approval, the user assumes all responsibility and liability for biological and chemical decontamination and cleansing. Labconco reserves the right to refuse delivery of any products which do not appear to have been properly cleaned and/or decontaminated prior to return.

Accessories	Part #	Description
	78169	Lower Deck Hexagonal Rotor Capacity for (18) 1.5ml microcentrifuge tubes 3ml Rotor (12mm x 75mm) Capacity for (72) 1.5ml centrifuge tubes or (16) 12 x 55 mm tubes or (32) 12 x 75 mm tubes
	78142	Rotor Capacity for (24) 1.5 ml microcentrifuge tubes or (6) 28 x 135 mm tubes or (6) 28 x 140 mm tubes
	78377	Rotor Capacity for (48) 1.5 microcentrifuge tubes or (16) 12 x 55 mm tubes or (16) 12 x 75 mm tubes or (12) 16 x 100 mm tubes or (12) 16 x 120 mm tubes
	78156	Glass Insert For Cold Trap
	78148	Acid Trap Insert
	78149	Moisture Trap Insert
	78150	Radioisotope Trap Insert
	78152	Solvent Trap Insert
	14721	Direct Drive Vacuum Pump 113 liter/minute pumping capacity with gas ballast. Blank-off pressure 2×10^{-4} Torr (0.2 micron). 115 VAC, 60 Hz, single phase, 6.3 amp

Accessory Part #	Description
14677	Direct Drive Vacuum Pump 190 Liter/minute pumping capacity with gas ballast. Blank off pressure 2×10^{-4} Torr (0.2 micron) (115 VAC, 60 Hz, single phase, 6.3 amp
77394	Direct Drive Vacuum Pump 113 liter/minute. Same as 14721 except 220/208-230 VAc, 50/60 Hz, single phase, 6.3 amp
77394-01	Direct Drive Vacuum Pump 190 liter/minute. Same as 14677 except 220/208-230 VAC, 50/60 Hz, single phase, 3.4 amp operation
79182	Vacuum Pump Diaphragm pump - corrosion resistant 115 V, 60 Hz 6.2 amps, single phase, 88 liter/min, maximum vacuum 29 in. Hg.
78182-01	Vacuum Pump Diaphragm pump - corrosion resistant 230 V, 50 Hz, 2.2 amps, single phase, 74 liter/min, maximum vacuum 29 in. Hg.
14723	Pump Exhaust Filter For direct drive vacuum pumps, not diaphragm pumps

If you have any questions that are not addressed in this manual, or if you need technical assistance, please contact Labconco's Customer Service Department at either (800) 821-5525, or (816) 333-8811, between the hours of 8:00 a.m. and 5:00 p.m. Central Standard Time.