


# Zero Pressure Weathercap, 12" Diameter

Model: # 5622100

## Description:

- PVC cap adds little static pressure to the exhaust system.
- Allows for vertical discharge of effluent air for dispersion away from the building.
- Equivalent resistance of straight duct: 5



Side View

Zero Pressure  
Weathercap for:



**LABCONCO CORPORATION** • KANSAS CITY MO USA 64132  
800-821-5525 / 816-333-8811 • labconco.com

Date:  
01-24-2025

Drawing #  
5622100

Drawn By:  
C. Alumbaugh