

# WaterPro® BT™ Systems

Buy a WaterPro® BT™ System and receive a **FREE** Filter Pack worth over \$1,040



Purchase a **WaterPro BT Purification System** now through December 31, 2026 and receive one Filter Pack **FREE**. Please see redemption form (on page 2).

9019200 WaterPro BT Filter Pack required for operation.


Filter Pack includes twin cartridges with pretreatment, reverse osmosis, activated carbon and deionization filtration. Shipping weight 6.6 lbs. (3 kg)



*WaterPro BT Systems produce both RO-purified and Type I water packaged in a compact all-in-one benchtop design.*

## Features:

- Ideal for laboratories needing 1-10 liters per day
- Delivers ultrapure Type 1 water of up to 18.2 megohm-cm at a typical rate of >0.5 liter per minute
- Small footprint design for benchtop or optional wall mounting
- Optional UV lamp further reduces bacteria and TOC levels necessary for life science applications
- Lightweight and easily installed or moved to share between labs



**LABCONCO**

# WaterPro® BT™ Filter Pack Redemption Form

## QUALIFYING PRODUCTS

Electrical	Description	Labconco Catalog Number
115V	WaterPro BT	9015020
230V North America	WaterPro BT	9015040
115V	WaterPro BT with Remote Dispenser	9015220
230V North America	WaterPro BT with Remote Dispenser	9015240
115V	WaterPro BT with UV light	9015120
230V North America	WaterPro BT with UV light	9015140
115V	WaterPro BT with UV and Remote Dispenser	9015320
230V North America	WaterPro BT with UV and Remote Dispenser	9015340

To redeem your free Filter Pack, complete the form below and fax to **(816) 363-0130** or e-mail to [promotions@labconco.com](mailto:promotions@labconco.com), attention: Customer Service.

\* indicates required fields

Company\* \_\_\_\_\_ Department \_\_\_\_\_  
First Name\* \_\_\_\_\_ Last Name\* \_\_\_\_\_  
Address\* \_\_\_\_\_  
Address 2 \_\_\_\_\_ Daytime Phone\* \_\_\_\_\_  
City\* \_\_\_\_\_ State\* \_\_\_\_\_ Zip Code\* \_\_\_\_\_  
Email\* \_\_\_\_\_ Fax \_\_\_\_\_  
Purchase Order # used\* \_\_\_\_\_ Date of Order\* \_\_\_\_\_  
Purchased from (name of distributor)\* \_\_\_\_\_

**Terms and Conditions:** Valid in US & Canada and select International Markets. One (1) redemption form must be completed for each WaterPro BT purchased. No limit of number of products purchased. Offer is valid on purchases made from January 1, 2026 through December 31, 2026. Redemption forms must be received by January 31, 2027. This offer is void where restricted or prohibited by law or company policy. May not be valid in combination with other Labconco offers unless stated in writing by a Labconco representative. Purchases made in countries other than the U.S. or Canada are not eligible. Regular shipping charges are included. Product appearance and specifications are subject to change without notice. Allow 4 to 8 weeks after receipt of completed redemption form for delivery of free Filter Pack.

### LABCONCO CORPORATION

8811 Prospect Avenue  
Kansas City, MO 64132  
(800) 821-5525 | (816) 363-0130 fax  
[labconco.com](http://labconco.com)

**LABCONCO**